ENGLISH

BODY MASSAGER

The massager vibration system affects effectively on muscular tissue, removes excessive fat and helps to get rid of cellulite.

The massager not only performs anti-cellulite function but also reduces backache and neck ache, helps to improve blood circulation and generally tone up the body, that's why the massager is perfectly suitable for people who spend a lot of time in sitting or standing position: office employees, teachers, doctors, salespersons and other staff.

Regular massage stimulates nerve endings, general strengthens your health, calms and helps to relieve stress.

DESCRIPTION

- Motor unit
- Massage intensity control knob
- Massage head
- Attachment for acupuncture massage
- Attachment for intensive massage
- Attachment for delicate massage
- Rolling attachment
- Protective attachment

SAFETY MEASURES

Before using the unit, read these instructions carefully. Keep these instructions for future reference. Use the unit according to its intended purpose only, as it is stated in this user manual. Mishandling the unit can lead to its breakage and cause harm to the user or damage to his/hers property.

Before switching the unit on, make sure that the operating voltage of the unit corresponds to voltage of your home mains.

- Continuous operation time while using the massager should not exceed 15 minutes; use the unit at least a 15-minute break between massage sessions.
- Do not use the unit near heat sources or open fire.
- Make sure that the power cord does not touch hot surfaces or sharp edges of furniture.
- Do not use the massager outdoors.
- Do not use the unit when you are sick or drowsy.
- Massage should give only pleasant sensations. If you feel pain or discomfort during massage, stop using the massager immediately and consult your doctor.
- Do not leave the operating unit unattended.
- Always switch the massager off and unplug it before cleaning and when you do not use it.
- Do not allow children to use the unit as a toy. For children safety reasons do not leave polyethylene bags, used as a packaging, unattended.
- Attention!** Do not allow children to play with polyethylene bags or packaging film. Danger of suffocation!
- This unit is not intended for usage by children or disabled persons unless they are given all the necessary and understandable instructions by a person who is responsible for their safety on safety measures and information about danger that can be caused by improper usage of the unit.
- Do not use the massager if your skin is damaged at the massage area.
- Do not use the unit, if the power cord, the plug or the unit body is damaged.
- Do not repair the unit by yourself. If the unit is damaged, apply to an authorized service center.

ATTENTION!

Before using the massager consult your doctor, if you have:

- cardiovascular diseases;
- the installed pacemaker;
- oncological diseases;
- spine implants or surgery;
- active chronic diseases;
- damaged skin or a skin disease in the massage area;

and if you:

- are pregnant;
- have recently undergone a surgery;
- are under treatment;
- have physical activity limits;

INDICATIONS FOR MASSAGER USE

You can use the massager for massaging practically every part of your body:

- shoulders (relieves tension in the neck or spine);
- forearms (tones up muscles);
- abdominals (tones up abdominal muscles and shapes the waistline);
- hips (tones up glutes);
- upper thighs (relieves tension in muscles of thigh and makes your legs slim);
- calf muscles (helps to take off weariness and improves blood circulation in legs blood vessels).

MASSAGE ATTACHMENTS (4, 5, 6, 7)

- Attachment for acupuncture massage (4) provides point massage and deep tissue massage.
- Attachment for intensive massage (5) improves blood circulation in muscles.
- Attachment for delicate massage (6) can be used for oil massage.
- Rolling attachment (7) is intended for massaging neck and shoulders.

Attention! Do not use the protective attachment (8) when using the rolling attachment (7).

Protective attachment (8)

Protective attachment (8) is intended for usage in areas where hair may get caught between rotating parts of the massage head (when massaging neck or shoulders).

Attention!

- Protective attachment (8) is not intended for usage with rolling attachment (7) and for oil massage.
- When using the massager without protective attachment (8) make sure that your hair does not get caught between rotating parts of the massage head (3).

BEFORE THE FIRST USE

After transportation or storage at low temperature it is necessary to keep it for at least two hours at room temperature before switching on.

- Unpack the unit and remove any stickers that can prevent unit operation.
- Check the unit for integrity; do not use the unit in case of damages.
- Before switching the unit on make sure that your home mains voltage corresponds to unit operating voltage.
- Wipe the motor unit (1) with a soft dry cloth.
- Wash the massage attachments (4, 5, 6, 7) and the protective attachment (8) with warm water and a neutral detergent and dry them thoroughly.

MASSAGING

- Attention!** Avoid direct impact on the vertebrae and perform all massage techniques on both sides of the spine during spine massage.
- Set the required massage attachment (4..7).
- Plug the power plug (3) into the mains socket.
- If you use the attachments (4, 5, 6), install the protective attachment (8) on the motor unit (1).

Note: Protective attachment (8) is not intended for usage with rolling attachment (7) and for oil massage.

- Insert the power plug into the mains socket.
- Before connecting the unit to the mains make sure that the massage intensity control knob (2) is set to the position «OFF».
- Turn the knob (2) to switch the massager on and select the required massage intensity.
- Massage the selected body area no longer than 3 minutes, then move the massager to another body area.
- After 15 minutes of massage, switch the unit off by setting the knob (2) to the position «OFF», and unplug it.
- Make at least a 15-minute break after every massage session.

RECOMMENDATIONS

- Use the massager in at least an hour after meals.
- Do not use the unit for head massage.
- Do not use the massager at ambient temperature above 35°С and in rooms with high humidity.
- Do not use the massager if you feel unwell.
- If you feel unpleasant symptoms during the massage, stop using the unit immediately and see your doctor.
- For massaging put the massager to the wanted body area with no excessive effort.
- Move the massager from the outside to the center of the body area.
- It is recommended to start massaging from cervical spine. Then you can proceed to massage of the back. Recommended massage time for one body area is from 3 to 5 minutes.
- You can use the massager for preventive anticellulite massage of buttocks or hips. Recommended duration of anticellulite massage session is from 4 to 5 minutes.
- Massage the abdominal area only clockwise.
- Do not repeat the unit by yourself. If the unit is damaged, apply to an authorized service center.
- Finish the session with feet massage; it activates acupuncture points and effects favorably on human body.
- General massage session should be no longer than 15 minutes with a subsequent 15-minute break.

CLEANING AND CARE

- Switch the unit off by setting the massage intensity control knob (2) to the position «OFF» and unplug the unit.
- Wipe the motor unit (1) with a soft dry cloth; if the motor unit (1) body is very dirty, wipe it with a slightly damp cloth with neutral detergent and then wipe it dry.
- Do not use abrasive substances and solvents to clean the motor unit (1).
- Do not immerse the motor unit (1), the power cord and the plug into water or any other liquids.
- Keep the unit away from children in a dry cool place.

DELIVERY SET

- Motor unit – 1 pc.
- Massage attachments – 4 pcs.
- Protective attachment (8) – 1 pc.
- Instruction manual – 1 pc.

TECHNICAL SPECIFICATIONS

Power supply: 220-240 V – 50 Hz
Power consumption: 50 W
Continuous operation time: not longer than 15 minutes

The manufacturer preserves the right to change the specifications of the unit without a preliminary notification.

Unit operating life is 5 years

Guarantee

Details regarding guarantee conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

This product conforms to the EMC-Requirements as laid down by the Council Directive 2004/108/EC, and to the Low Voltage Regulation (2006/95/EC)

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