Contents

English 01-04
Deutsch 05-08
Français 09-12
Español 13-16
Italiano 17-20

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Introduction

Thank you for purchasing the Nairo Foot Massager with Kneading Massage and Heat. Please read this user manual carefully and keep it in a safe place for future reference. If you need any assistance, please contact our support team with your product model number and Amazon order number.

Included

1 x Foot massager
1 x User manual
1 x Adapter

Product Specifications

Model	MGF-1005
Name	Nairo Foot Massager with Kneading Massage and Heat
Voltage	12V 1.5A
Net weight	3.5lbs / 1.6kg
Default timer	15:11 minutes
Functions	Kneading / Vibrating / Heating

Product Diagram

Important Safeguards

- Please do not use the product in a humid environment (like a bathroom) to avoid electric shocks or other safety problems.
- Please do not expose the product to direct sunlight or high temperature.
- Please do not sit on the product.
- Please avoid the product being damaged by sharp objects. Do not heat the power cable. Do not use the product if the power cable is damaged.
- Please disconnect the product from the power supply before cleaning. Do not clean with corrosive detergents.
- Please store the product in a safe, dry, and airy place; if the product is not to be used for long, please cover it with dust-proof cloth to prevent the accumulation of dust on it.

Q & A

Problem	Cause	Solution
There exists abnormal smell.	The product will reach the customer directly after leaving the factory. A smell of plastic will remain after being packed for long.	Put the product in an airy place without being packed for several days to eliminate the abnormal smell.
The product works noisily.	The product will cause noise when working with vibrating function.	The noise could be somewhat abated by placing the product on a flat surface.
The product stops working automatically.	The product will shut off automatically after 15 minutes of constant use.	Press the "ON/OFF" button to restart the product.

Product Operation Guide

Kneading massage, heat, and vibration are available. After connected to the power supply, the product will shift to the standby mode. Below is the instruction of its features:

1. The "ON/OFF" Button

After connected to the power supply, the product will shift to the standby mode. Press this button to turn on the product with the corresponding indicator on and it starts being operational. The kneading massage mode I, the heating function and the heating function would be activated by default with the timer beginning to work (15 minutes). When the product is on, press this button, and the product will shift back to the standby mode. The product will also shift back to the standby mode if time is up.

2. The "Mode" Button (2 massage directions / automatic direction changing is optional)

If the product works without the kneading massage, press this button to activate kneading massage mode I with the mode indicator remaining on all along. When the product is working with kneading massage mode I, press this button to shift to kneading massage mode II with the mode indicator flickering slowly. When the product is working with kneading massage mode II, press this button to shift to kneading massage mode III with the mode indicator flickering rapidly. When the product is working with kneading massage mode III, press this button to stop the kneading massage with the mode indicator off.

3. The "Vibration" Button

If the product works without the vibrating function, press this button to activate the vibrating function with the corresponding indicator on; if the vibrating function is on, press this button to turn off the vibrating function with the corresponding indicator off.

4. The "Heat" Button

If the product works without the heating function, press this button to activate the heating function with the corresponding indicator on; if the heating function is on, press this button to turn off the heating function with the corresponding indicator off.

Warranty & Customer Support

For questions, assistance, support, or warranty claims, please contact us via the email or the address below:
Web: www.naipocare.com
Email: support@naipocare.com
Address: Huanan City Electronic Trading Building P07 Room 120, Longgang District, Shenzhen, China.

*Note: Products purchased directly from Nairo benefit from our out-of-idea service. If the purchases are made from other selling teams, we recommend you contact them for warranty issues.

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